



PRAF1 rabbit pAb

Catalog No	YP-Ab-09696
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	RABAC1 PRA1 PRAF1
Immunogen	Synthesized peptide derived from human PRAF1 AA range: 130-180
Specificity	This antibody detects endogenous levels of human PRAF1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	A49;DNA directed RNA polymerase I subunit RPA49;DNA-directed RNA polymerase I subunit E;DNA-directed RNA polymerase I subunit RPA49;FLJ 1339;FLJ1397;PAF53;POLR 1 E;POLR1E;Polymerase (RNA) I associated factor 1;Polymerase (RNA) I polypeptide E;53kDa;PRAF1;RNA polymerase I associated factor 1;RNA polymerase I associated factor 53;RNA polymerase I polypeptide E;RNA polymerase I subunit A49;RNA polymerase I-associated factor 1;RNA polymerase I-associated factor 53;RP1145L18.3;RPA49_HUMAN; RABAC1 PRA1 PRAF1;Prenylated Rab acceptor protein 1 (PRA1 family protein 1)
Observed Band	
Calculated Molecular Weight	54 kDa; 47kD
Cell Pathway	Cell membrane ; Multi-pass membrane protein . Cytoplasm . Golgi apparatus . Cytoplasmic vesicle, secretory vesicle, synaptic vesicle . According to some authors, it is an integral membrane protein, while others showed that it is cytoplasmic and membrane-associated to Golgi and synaptic vesicles. .
Tissue Specificity	Ubiquitous. Strongest expression found in placenta, pituitary gland, kidney, lung and stomach.
Function	



Background

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images