



FSCN2 rabbit pAb

Catalog No	YP-Ab-09697
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	FSCN2
Immunogen	Synthesized peptide derived from human FSCN2 AA range: 356-406
Specificity	This antibody detects endogenous levels of human FSCN2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Fascin homolog 2;actin bundling protein;retinal Strongylocentrotus purpuratus; Fascin Strongylocentrotus purpuratus homolog 2 actin bundling protein;retinal; Fascin-2;FSCN2;FSCN2_HUMAN;Retinal fascin;RFSN;RP3;Fascin-2 (Retinal fascin)
Observed Band	
Calculated Molecular Weight	55kD
Cell Pathway	Cytoplasm, cytoskeleton . Cell projection, stereocilium .
Tissue Specificity	Localized specifically in the outer and inner segments of the photoreceptor cells in the retina.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
