



FHL17 rabbit pAb

Catalog No	YP-Ab-09698
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	FTHL17
Immunogen	Synthesized peptide derived from human FHL17 AA range: 67-117
Specificity	This antibody detects endogenous levels of human FHL17
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FHL17;FTHL17;Ferritin heavy polypeptide-like 17 (Cancer/testis antigen 38) (CT38)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	Testis specific. Also expressed in several cancers.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
