



ELP6 rabbit pAb

Catalog No	YP-Ab-09699
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ELP6 ATP1 C3orf75 TMEM103
Immunogen	Synthesized peptide derived from human ELP6 AA range: 22-72
Specificity	This antibody detects endogenous levels of human ELP6
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ELP6;ELP6 ATP1 C3orf75 TMEM103;Elongator complex protein 6 (Angiotonin-transactivated protein 1) (Protein TMEM103)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
