



CUTC rabbit pAb

Catalog No	YP-Ab-09702
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	CUTC CGI-32
Immunogen	Synthesized peptide derived from human CUTC AA range: 174-224
Specificity	This antibody detects endogenous levels of human CUTC
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CGI-32;CUTC;CUTC CGI-32;Copper homeostasis protein cutC homolog
Observed Band	29 kDa
Calculated Molecular Weight	29 kDa
Cell Pathway	Cytoplasm . Nucleus . The overexpressed protein is detected in the cytoplasm, and depending on the cell line, also in the nucleus.
Tissue Specificity	Ubiquitous.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
