



STRBP rabbit pAb

Catalog No	YP-Ab-09703
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	STRBP SPNR
Immunogen	Synthesized peptide derived from human STRBP AA range: 255-305
Specificity	This antibody detects endogenous levels of human STRBP
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SPNR;Spermatid perinuclear RNA-binding protein;p74;ILF3L;STRBP;STRBP SPNR
Observed Band	74 kDa
Calculated Molecular Weight	672 aa, 74 kDa
Cell Pathway	Cytoplasm . Microtubule-associated that localizes to the manchette in developing spermatids. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
