



CYTSB rabbit pAb

Catalog No	YP-Ab-09709
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	SPECC1 CYTSB NSP5
Immunogen	Synthesized peptide derived from human CYTSB AA range: 403-453
Specificity	This antibody detects endogenous levels of human CYTSB
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	cytospin B;CYTSB;HCMOGT 1;NSP;NSP5;Nuclear structure protein 5;SPECC1;Sperm antigen HCMOGT 1;SPECC1 CYTSB NSP5;Cytospin-B (Nuclear structure protein 5) (NSP5) (Sperm antigen HCMOGT-1) (Sperm antigen with calponin homology and coiled-coil domains 1)
Observed Band	119 kDa
Calculated Molecular Weight	119 kDa
Cell Pathway	Nucleus .; [Isoform 3]: Membrane ; Lipid-anchor .; [Isoform 4]: Membrane ; Lipid-anchor .; [Isoform 5]: Membrane ; Lipid-anchor .
Tissue Specificity	Highly expressed in testis. Barely detectable in other tissues. Also highly expressed in some cancer cell lines.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
