



ODFP4 rabbit pAb

Catalog No	YP-Ab-09710
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ODF4 OPPO1
Immunogen	Synthesized peptide derived from human ODFP4 AA range: 119-169
Specificity	This antibody detects endogenous levels of human ODFP4
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ODFP4;ODF4 OPPO1;Outer dense fiber protein 4 (Outer dense fiber of sperm tails protein 4) (Testis-specific protein oppo 1) (hOPPO1)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Multi-pass membrane protein .
Tissue Specificity	Expressed in testis and sperm; especially localized to sperm tail (at protein level)
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
