



TEX36 rabbit pAb

Catalog No	YP-Ab-09713
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	TEX36 C10orf122
Immunogen	Synthesized peptide derived from human TEX36 AA range: 102-152
Specificity	This antibody detects endogenous levels of human TEX36
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TEX36;TEX36 C10orf122;Testis-expressed sequence 36 protein
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
