



EVI2A rabbit pAb

Catalog No	YP-Ab-09715
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	EVI2A EVDA EVI2
Immunogen	Synthesized peptide derived from human EVI2A AA range: 135-185
Specificity	This antibody detects endogenous levels of human EVI2A
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EVI-2A;EVI2;EVI 2A;EVDA;EVI2A;EVI2A EVDA EVI2;Protein EVI2A (Ecotropic viral integration site 2A protein homolog) (EVI-2A)
Observed Band	26 kDa
Calculated Molecular Weight	236 aa, 26 kDa
Cell Pathway	Membrane; Single-pass type I membrane protein.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
