



ZNF74 rabbit pAb

Catalog No	YP-Ab-09716
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	ZNF74 ZNF520
Immunogen	Synthesized peptide derived from human ZNF74 AA range: 577-627
Specificity	This antibody detects endogenous levels of human ZNF74
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Cos52;hZNF7;Zfp520;zinc finger protein 74;ZNF520;ZNF74;ZNF74 ZNF520; Zinc finger protein 74 (Zinc finger protein 520) (hZNF7)
Observed Band	64-72 kDa
Calculated Molecular Weight	644 aa, 72 kDa
Cell Pathway	Nucleus.
Tissue Specificity	Highly expressed in the fetal brain.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
