



SPG20 rabbit pAb

Catalog No	YP-Ab-09717
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	SPG20 KIAA0610 TAHCCP1
Immunogen	Synthesized peptide derived from human SPG20 AA range: 303-353
Specificity	This antibody detects endogenous levels of human SPG20
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SPG20;KIAA0610;SPART;Spartin;SPG 20;SPG20 KIAA0610 TAHCCP1;Spartin (Spastic paraplegia 20 protein) (Trans-activated by hepatitis C virus core protein 1)
Observed Band	75-84 kDa
Calculated Molecular Weight	72 kDa
Cell Pathway	Cytoplasm . Midbody . Transiently associated with endosomes (PubMed:19580544). Colocalized with IST1 to the ends of Flemming bodies during cytokinesis (PubMed:20719964). .
Tissue Specificity	Ubiquitously expressed, with highest levels of expression detected in adipose tissue.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images