



GALP rabbit pAb

Catalog No	YP-Ab-09719
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	GALP
Immunogen	Synthesized peptide derived from human GALP AA range: 12-62
Specificity	This antibody detects endogenous levels of human GALP
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	galanin like peptide;GALP;Galanin-like peptide
Observed Band	5 kDa
Calculated Molecular Weight	116 aa, 13 kDa
Cell Pathway	Secreted.
Tissue Specificity	Isoform 2 is found in ganglia of ganglioneuroma and ganglioneuroblastoma, as well as in differentiated tumor cells of neuroblastoma tissues. Not found in undifferentiated neuroblasts. Isoform 2 is found in the skin, in pericytes covering microvascular arterioles and venules on their abluminal surfaces. In larger vessels, isoform 2 is expressed in layers of smooth muscle cells. Isoform 2 is not detected in endothelial cells.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images