



TRH rabbit pAb

Catalog No	YP-Ab-09721
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TRH
Immunogen	Synthesized peptide derived from human TRH AA range: 133-183
Specificity	This antibody detects endogenous levels of human TRH
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Prothyroliberin;Pro-thyrotropin-releasing hormone;Protirelin;Thyrotropin-releasing hormone;TRF;TRH;Prothyroliberin [Cleaved into: Thyroliberin (Protirelin) (TSH-releasing factor) (Thyrotropin-releasing factor) (TRF) (Thyrotropin-releasing hormone) (TRH)]
Observed Band	27-33 kDa
Calculated Molecular Weight	242 aa, 27 kDa
Cell Pathway	Secreted.
Tissue Specificity	Hypothalamus. Expressed in the hair follicle epithelium (at protein level).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
