



TRI10 rabbit pAb

Catalog No	YP-Ab-09722
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	TRIM10 RFB30 RNF9
Immunogen	Synthesized peptide derived from human TRI10 AA range: 319-369
Specificity	This antibody detects endogenous levels of human TRI10
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TRI10;TRIM10 RFB30 RNF9;Tripartite motif-containing protein 10 (B30-RING finger protein) (RING finger protein 9)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
