



UD2B7 rabbit pAb

Catalog No	YP-Ab-09725
Isotype	IgG
Reactivity	Human;Rat
Applications	WB
Gene Name	UGT2B7 UGTB2B9
Immunogen	Synthesized peptide derived from human UD2B7 AA range: 467-517
Specificity	This antibody detects endogenous levels of human UD2B7
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	UDP glucuronosyltransferase 2 family;member B7;3 4 catechol estrogen specific; 3;4 catechol estrogen specific UDPGT;4-catechol estrogen-specific UDPGT; UD2B7_HUMAN;polypeptide B7;UDP glucuronosyltransferase 2 family polypeptide B7;member B9;UDP glucuronosyltransferase 2B7;UDP glucuronosyltransferase2 family;UDP glucuronyltransferase family 2 beta 7;UDP glycosyltransferase 2 family polypeptide B7;UDP-glucuronosyltransferase 2B7; UDP-glucuronosyltransferase 2B9;UDPGT;UDPGT 2B7;UDPGT 2B9; UDPGT2B7;UDPGTh 2;UDPGTh-2;UDPGTh2;UGT2B7;UGT2B9;UGTB2B9; uridine diphosphate glycosyltransferase 2 family;UD2B7;UGT2B7 UGTB2B9; UDP-glucuronosyltransferase 2B7 (UDPGT 2B7) (EC 2.4.1.17) (3;4-catechol estrogen-specific UDPGT) (UDP-glucuronosyltransferase 2B9) (UDPGT 2B9) (UDPGTh-2)
Observed Band	
Calculated Molecular Weight	55 kDa; 61kD
Cell Pathway	Endoplasmic reticulum membrane ; Single-pass membrane protein .
Tissue Specificity	
Function	



Background

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
