



DJB14 rabbit pAb

Catalog No	YP-Ab-09728
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	DNAJB14 UNQ9427/PRO34683
Immunogen	Synthesized peptide derived from human DJB14 AA range: 159-209
Specificity	This antibody detects endogenous levels of human DJB14
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	573496F1Rik;DJB14_HUMAN;DnaJ (Hsp4) homolog;subfamily B;member 14; DnaJ homolog subfamily B member 14;dnajb14;EGNR9427;FLJ14281; MGC22187;PRO34683;UNQ9427/PRO34683;DJB14;DNAJB14 UNQ9427/ PRO34683
Observed Band	
Calculated Molecular Weight	43 kDa; 43kD
Cell Pathway	Endoplasmic reticulum membrane ; Single-pass membrane protein . Nucleus membrane ; Single-pass membrane protein . Localizes to the endoplasmic reticulum membrane (PubMed:23018488, PubMed:24732912, PubMed:27916661). When overexpressed, forms membranous structures in the nucleus (PubMed:24732912). .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images