



XKR2 rabbit pAb

Catalog No	YP-Ab-09731
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	XKRX XKR2 XPLAC XRG2
Immunogen	Synthesized peptide derived from human XKR2 AA range: 376-426
Specificity	This antibody detects endogenous levels of human XKR2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Membrane protein XPLAC;X Kell blood group-related;X-linked;XK related protein 2;XKR2;XK-related protein 2;XKRX XKR2 XPLAC XRG2;XK-related protein 2 (Membrane protein XPLAC) (X Kell blood group-related;X-linked)
Observed Band	30 kDa
Calculated Molecular Weight	462 aa, 54 kDa
Cell Pathway	Cell membrane ; Multi-pass membrane protein .
Tissue Specificity	Expressed predominantly in the placenta, in syncytiotrophoblasts. Moderate levels in the adrenal gland, low levels in the trachea and very low levels in the bone marrow.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
