



COX7S rabbit pAb

Catalog No	YP-Ab-09732
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	COX7A2P2 COX7A3 COX7AL2 COX7AP2
Immunogen	Synthesized peptide derived from human COX7S AA range: 13-63
Specificity	This antibody detects endogenous levels of human COX7S
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	COX7A3;COX7AL2;COX7S/A2;Cytochrome c oxidase subunit 7A3; mitochondrial;Cytochrome c oxidase subunit VIIa 3;COX7A2;COX7AL; COX7AL1;COXVIIa L;COXVIIAL;CX7A2_HUMAN;Cytochrome c oxidase polypeptide 7A2 mitochondrial;Cytochrome c oxidase polypeptide VIIa liver/ heart;Cytochrome c oxidase subunit 7A2;Cytochrome c oxidase subunit 7A2 mitochondrial;Cytochrome c oxidase subunit VIIa polypeptide 2;Cytochrome c oxidase subunit VIIa-L;Cytochrome c oxidase subunit VIIa-liver/heart; Cytochrome c oxidase subunit VIIaL;Hepatic cytochrome c oxidase chain VIIa; MGC11895;MGC118951;MGC118952;MGC126875;MGC126877;VIIAL;COX7S; COX7A2P2 COX7A3 COX7AL2 COX7AP2;Putative cytochrome c oxidase subunit 7A3;mitochondrial (Cytochrome c oxidase subunit VIIa 3)
Observed Band	
Calculated Molecular Weight	12kDa; 12kD,9kD
Cell Pathway	Mitochondrion inner membrane .
Tissue Specificity	
Function	



Background

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
