



TIM44 rabbit pAb

Catalog No	YP-Ab-09734
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TIMM44 MIMT44 TIM44
Immunogen	Synthesized peptide derived from human TIM44 AA range: 260-310
Specificity	This antibody detects endogenous levels of human TIM44
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TIMM 44;TIM44;MIMT44;TIMM44 MIMT44 TIM44;Mitochondrial import inner membrane translocase subunit TIM44
Observed Band	45 kDa
Calculated Molecular Weight	452 aa, 51 kDa
Cell Pathway	Mitochondrion inner membrane ; Peripheral membrane protein ; Matrix side . Mitochondrion matrix .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
