



TR16L rabbit pAb

Catalog No	YP-Ab-09735
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	TRIM16L TRIM70
Immunogen	Synthesized peptide derived from human TR16L AA range: 91-141
Specificity	This antibody detects endogenous levels of human TR16L
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TR16L;TRIM16L TRIM70;Tripartite motif-containing protein 16-like protein (Tripartite motif-containing protein 70)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
