



CF089 rabbit pAb

Catalog No	YP-Ab-09736
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	C6orf89 BRAP UNQ177/PRO203
Immunogen	Synthesized peptide derived from human CF089 AA range: 193-243
Specificity	This antibody detects endogenous levels of human CF089
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CF089;C6orf89 BRAP UNQ177/PRO203;Bombesin receptor-activated protein C6orf89 (Amfion)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	[Isoform 1]: Golgi apparatus membrane ; Single-pass type II membrane protein . Midbody . ; [Isoform 2]: Golgi apparatus membrane ; Single-pass type II membrane protein . Midbody . ; [Isoform 3]: Cytoplasm . Nucleus, nucleolus . Retained in nucleolar organiser regions (NORs) in mitotic cells (PubMed:23460338). .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images