



BTNL3 rabbit pAb

Catalog No	YP-Ab-09737
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	BTNL3 BTNLR COLF4100 UNQ744/PRO1472
Immunogen	Synthesized peptide derived from human BTNL3 AA range: 269-319
Specificity	This antibody detects endogenous levels of human BTNL3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	BTNL3;BTNL3 BTNLR COLF4100 UNQ744/PRO1472;Butyrophilin-like protein 3 (Butyrophilin-like receptor)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Single-pass type I membrane protein .
Tissue Specificity	Expressed in small intestine, colon, testis, spleen, and leukocyte.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
