



DDAH2 rabbit pAb

Catalog No	YP-Ab-09739
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	DDAH2 DDAH G6A NG30
Immunogen	Synthesized peptide derived from human DDAH2 AA range: 179-229
Specificity	This antibody detects endogenous levels of human DDAH2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DDAH;DDAH 2;DDAH-2;DDAHII;Dimethylargininase 2;DDAH2;DDAH2 DDAH G6A NG30;N(G);N(G)-dimethylarginine dimethylaminohydrolase 2 (DDAH-2) (Dimethylarginine dimethylaminohydrolase 2) (EC 3.5.3.18) (DDAHII) (Dimethylargininase-2) (Protein G6a) (S-phase protein)
Observed Band	30 kDa
Calculated Molecular Weight	29 kDa
Cell Pathway	Cytoplasm . Mitochondrion . Translocates from cytosol to mitochondrion upon IL-1beta stimulation in chondrocytes.
Tissue Specificity	Detected in heart, placenta, lung, liver, skeletal muscle, kidney and pancreas, and at very low levels in brain.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images