



PALD rabbit pAb

Catalog No	YP-Ab-09744
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	PALD1 KIAA1274 PALD
Immunogen	Synthesized peptide derived from human PALD AA range: 266-316
Specificity	This antibody detects endogenous levels of human PALD
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PALD;PALD1 KIAA1274 PALD;Paladin
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm, cytosol .
Tissue Specificity	Expressed in endothelial cells, and in certain larger vessels, in mural cells. In the brain, possibly expressed in microglia. Expressed in peripheral blood mononuclear cells (at protein level).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

