



NCKPL rabbit pAb

Catalog No	YP-Ab-09762
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	NCKAP1L HEM1
Immunogen	Synthesized peptide derived from human NCKPL AA range: 582-632
Specificity	This antibody detects endogenous levels of human NCKPL
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	493568P13Rik;AI46383;HEM1;Hematopoietic protein 1;HEMATOPOIETIC PROTEIN HEM-1;Hemp1;Membrane associated protein hem1;Membrane-associated protein HEM-1;NCK associated protein 1 like;Nck-associated protein 1-like;NCKAP1L;NCKPL_HUMAN;NCKPL;NCKAP1L HEM1;Nck-associated protein 1-like (Hematopoietic protein 1) (Membrane-associated protein HEM-1)
Observed Band	124kD
Calculated Molecular Weight	124kD
Cell Pathway	Cell membrane ; Single-pass membrane protein ; Cytoplasmic side . Cytoplasm . Localizes to the leading edge of polarized neutrophils. .
Tissue Specificity	Expressed only in cells of hematopoietic origin (PubMed:7643388, PubMed:1932118). Expressed in neutrophils (at protein level) (PubMed:16417406). Expressed in T-cells (at protein level) (PubMed:32647003).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images