



ZNF8 rabbit pAb

Catalog No	YP-Ab-09763
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	ZNF8
Immunogen	Synthesized peptide derived from human ZNF8 AA range: 27-77
Specificity	This antibody detects endogenous levels of human ZNF8
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Zfp128;Zinc finger protein 8;ZNF8;Zinc finger protein 8 (Zinc finger protein HF.18)
Observed Band	65 kDa
Calculated Molecular Weight	575 aa, 65 kDa
Cell Pathway	Nucleus .
Tissue Specificity	Ubiquitously present in many human cell lines of different embryological derivation.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
