



CL18C rabbit pAb

Catalog No	YP-Ab-09765
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	CLEC18C MRLP3
Immunogen	Synthesized peptide derived from human CL18C AA range: 272-322
Specificity	This antibody detects endogenous levels of human CL18C
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CL18C;CLEC18C MRLP3;C-type lectin domain family 18 member C (Mannose receptor-like protein 3)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Secreted . Endoplasmic reticulum . Golgi apparatus . Endosome .
Tissue Specificity	Detected in peripheral blood cells.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
