



COFA1 rabbit pAb

Catalog No	YP-Ab-09767
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	COL15A1
Immunogen	Synthesized peptide derived from human COFA1 AA range: 1231-1280
Specificity	This antibody detects endogenous levels of human COFA1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	COFA1;COL15A1;Collagen alpha-1(XV) chain [Cleaved into: Endostatin (Endostatin-XV) (Related to endostatin) (Restin)]
Calculated Molecular Weight	153kD
Cell Pathway	Secreted, extracellular space, extracellular matrix .
Tissue Specificity	Expressed predominantly in internal organs such as adrenal gland, pancreas and kidney.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

