



COX14 rabbit pAb

Catalog No	YP-Ab-09782
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	COX14 C12orf62
Immunogen	Synthesized peptide derived from human COX14 AA range: 4-54
Specificity	This antibody detects endogenous levels of human COX14
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C12orf62;COX14;COX14 C12orf62;Cytochrome c oxidase assembly protein COX14
Observed Band	7-12 kDa
Calculated Molecular Weight	
Cell Pathway	Mitochondrion membrane ; Single-pass membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
