



NETR rabbit pAb

Catalog No	YP-Ab-09786
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PRSS12
Immunogen	Synthesized peptide derived from human NETR AA range: 275-325
Specificity	This antibody detects endogenous levels of human NETR
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	epididymis tissue protein Li 131P;NET1_HUMAN;Netrin 1;mouse;homolog of; Netrin 1-like;Netrin-1;Netrin1;Ntn 1;NTN 1L;NTN1;NTN1L;Unc6;NETR;PRSS12; Neurotrypsin (EC 3.4.21.-) (Leydin) (Motopsin) (Serine protease 12)
Observed Band	96kD
Calculated Molecular Weight	67 kDa; 68kD
Cell Pathway	Secreted.
Tissue Specificity	Brain and Leydig cells of the testis.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
