



K2C80 rabbit pAb

Catalog No	YP-Ab-09793
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	KRT80 KB20
Immunogen	Synthesized peptide derived from human K2C80 AA range: 107-157
Specificity	This antibody detects endogenous levels of human K2C80
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CK 8;CK-8;CK8;Cytokeratin 8;Cytokeratin-8;K 8;K2C8_HUMAN;K8;KB2;Keratin 8;Keratin;Keratin-8;KRT 8;Krt8;type II cytoskeletal 8;Type II keratin;Type-II keratin Kb2;K2C80;KRT80 KB20;type II cytoskeletal 80 (Cytokeratin-80) (CK-80) (Keratin-80) (K80) (Type-II keratin Kb20)
Observed Band	
Calculated Molecular Weight	5kD; 51kD
Cell Pathway	
Tissue Specificity	Weakly expressed in tongue, but not skin or in any other tissues or organs examined.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

