



UN13D rabbit pAb

Catalog No	YP-Ab-09794
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	UNC13D
Immunogen	Synthesized peptide derived from human UN13D AA range: 213-263
Specificity	This antibody detects endogenous levels of human UN13D
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FHL 3;FHL3;FLJ67;HLH 3;HLH3;HPLH 3;HPLH3;Jinx;Munc13 4;Munc13-4; Protein unc 13 homolog D;Protein unc-13 homolog D;UN13D_HUMAN;Unc 13 homolog D;UNC 13D;Unc-13 homolog D (C. elegans);Unc13 homolog D (C elegans);Unc13 homolog D;UNC13;C. elegans;homolog of;D;UNC13D;Unc13h4; UN13D;Protein unc-13 homolog D (Munc13-4)
Observed Band	
Calculated Molecular Weight	123 kDa; 123kD
Cell Pathway	Cytoplasm. Membrane; Peripheral membrane protein. Late endosome. Recycling endosome. Lysosome. Colocalizes with cytotoxic granules at the plasma membrane. Localizes to endosomal exocytic vesicles.
Tissue Specificity	Expressed at high levels in spleen, thymus and leukocytes. Also expressed in lung and placenta, and at very low levels in brain, heart, skeletal muscle and kidney. Expressed in cytotoxic T-lymphocytes (CTL) and mast cells.
Function	
Background	

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images