



HEMH rabbit pAb

Catalog No	YP-Ab-09804
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	FECH
Immunogen	Synthesized peptide derived from human HEMH AA range: 169-219
Specificity	This antibody detects endogenous levels of human HEMH
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EPP;FCE;Fch;Fcl;fech;Ferrochelatase (protoporphyrin);Ferrochelatase; Ferrochelatase mitochondrial;Heme synthase;Heme synthetase;HEMH_HUMAN; mitochondrial;Protoheme ferro lyase;Protoheme ferro-lyase;Protoporphyrin; HEMH;mitochondrial (EC 4.99.1.1) (Heme synthase) (Protoheme ferro-lyase)
Observed Band	
Calculated Molecular Weight	46 kDa; 48kD
Cell Pathway	Mitochondrion inner membrane; Peripheral membrane protein; Matrix side.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
