



CB081 rabbit pAb

Catalog No	YP-Ab-09809
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	C2orf81
Immunogen	Synthesized peptide derived from human CB081 AA range: 433-483
Specificity	This antibody detects endogenous levels of human CB081
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CB081;C2orf81;Uncharacterized protein C2orf81
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
