



ERI3 rabbit pAb

Catalog No	YP-Ab-09811
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	ERI3 PINT1 PRNPIP PRNPIP1
Immunogen	Synthesized peptide derived from human ERI3 AA range: 64-114
Specificity	This antibody detects endogenous levels of human ERI3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	3' 5' exonuclease ERI 1;3' 5' exonuclease ERI1;3' EXO;3' exoribonuclease;3' HEXO;3' histone mRNA exonuclease 1;3"-5" exonuclease eri1;3"-5" exoribonuclease 1;3'EXO;3'HEXO;ERI 1;Eri 1 homolog;Eri-1 homolog;eri1;Eri1 homolog;ERI1_HUMAN;HEXO;Histone mRNA 3' end specific exonuclease; Histone mRNA 3' end specific exoribonuclease;Histone mRNA 3"-end-specific exoribonuclease;Histone mRNA 3"-exonuclease 1;MGC35395;Protein 3"hExo; THEX 1;THEX1;Three prime histone mRNA exonuclease 1;ERI3;ERI3 PINT1 PRNPIP PRNPIP1;ERI1 exoribonuclease 3 (EC 3.1.-.-) (Prion interactor 1) (Prion protein-interacting protein)
Observed Band	
Calculated Molecular Weight	4 KD; 4kD
Cell Pathway	
Tissue Specificity	
Function	
Background	

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
