



PITM3 rabbit pAb

Catalog No	YP-Ab-09812
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	PITPNM3 NIR1
Immunogen	Synthesized peptide derived from human PITM3 AA range: 563-613
Specificity	This antibody detects endogenous levels of human PITM3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PITM3;PITPNM3 NIR1;Membrane-associated phosphatidylinositol transfer protein 3 (Phosphatidylinositol transfer protein;membrane-associated 3) (PITPnm 3) (Pyk2 N-terminal domain-interacting receptor 1) (NIR-1)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Endomembrane system ; Peripheral membrane protein .
Tissue Specificity	Detected in brain and spleen, and at low levels in ovary.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
