



WFKN1 rabbit pAb

Catalog No	YP-Ab-09813
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	WFIKKN1 C16orf12 GASP2 WFIKKN PSEC0040
Immunogen	Synthesized peptide derived from human WFKN1 AA range: 83-133
Specificity	This antibody detects endogenous levels of human WFKN1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	WFKN1;WFIKKN1 C16orf12 GASP2 WFIKKN PSEC0040;WAP;kazal; immunoglobulin;kunitz and NTR domain-containing protein 1 (Growth and differentiation factor-associated serum protein 2) (GASP-2) (hGASP-2) (WAP; follistatin;kunitz and N
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Secreted .
Tissue Specificity	Expressed in pancreas, kidney, liver, placenta, and lung.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
