



OR9A1 rabbit pAb

Catalog No	YP-Ab-09815
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	OR9A1P OR9A1
Immunogen	Synthesized peptide derived from human OR9A1 AA range: 65-115
Specificity	This antibody detects endogenous levels of human OR9A1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	OR9A1;OR9A1P OR9A1;Putative olfactory receptor 9A1 (HSHTPRX06)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell membrane; Multi-pass membrane protein.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

