



CTG1B rabbit pAb

Catalog No	YP-Ab-09826
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	CTAG1A CTAG CTAG1 ESO1 LAGE2 LAGE2A; CTAG1B LAGE2B
Immunogen	Synthesized peptide derived from human CTG1B AA range: 129-179
Specificity	This antibody detects endogenous levels of human CTG1B
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CTAG1A;CTAG1B;Autoimmunogenic cancer/testis antigen NY ESO 1; Autoimmunogenic cancer/testis antigen NY-ESO-1;Cancer antigen 3;Cancer/testis antigen 1;Cancer/testis antigen 1B;Cancer/testis antigen 6.1;CT6.1;CTAG 1;CTAG 1B;CTAG;CTAG1;CTG1B_HUMAN;ESO 1;ESO1;L antigen family member 2;LAGE 2;LAGE 2 protein;LAGE 2B;LAGE-2;LAGE2;LAGE2 protein;LAGE2A;LAGE2B;New York esophageal squamous cell carcinoma 1;NY ESO 1; NYESO 1;NYESO1;CTG1B;CTAG1A CTAG CTAG1 ESO1 LAGE2 LAGE2A; CTAG1B LAGE2B;Cancer/testis antigen 1 (Autoimmunogenic cancer/testis antigen NY-ESO-1) (Cancer/testis antigen 6.1) (CT6.1) (L antigen family member 2) (LAGE-2)
Observed Band	
Calculated Molecular Weight	18kD
Cell Pathway	Cytoplasm .
Tissue Specificity	Expressed in testis and ovary and in a wide variety of cancers. Detected in uterine myometrium. Expressed from 18 weeks until birth in human fetal testis. In the adult testis, is strongly expressed in spermatogonia and in primary spermatocytes, but not in post-meiotic cells or in testicular somatic cells (at protein level).

Function**Background****matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images