



ISM1 rabbit pAb

Catalog No	YP-Ab-09835
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	ISM1 C20orf82 ISM
Immunogen	Synthesized peptide derived from human ISM1 AA range: 312-362
Specificity	This antibody detects endogenous levels of human ISM1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	bA149I18.1;C20orf82;dJ1077I2.1;ISM;ISM1;Isthmin;Isthmin 1;isthmin 1 homolog (zebrafish);ISM1 C20orf82 ISM;Isthmin-1
Observed Band	52 kDa, 70 kDa
Calculated Molecular Weight	52kd
Cell Pathway	Secreted .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
