



QRIC2 rabbit pAb

Catalog No	YP-Ab-09842
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	QRICH2
Immunogen	Synthesized peptide derived from human QRIC2 AA range: 1251-1301
Specificity	This antibody detects endogenous levels of human QRIC2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	QRIC2;QRICH2;Glutamine-rich protein 2
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus membrane . Nucleus . Cytoplasm . Cell projection, cilium, flagellum . Localization varies during spermatozoa development. The protein is distributed in the nuclear membrane of the spermatogonia, in the nucleus of round spermatids, in the nucleus and cytoplasm of early elongating spermatids, in the cytoplasm of late elongating spermatids, and in the flagella of epididymal spermatozoa. .
Tissue Specificity	Expressed in the sperm.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images