



U2AFL rabbit pAb

Catalog No	YP-Ab-09854
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	ZRSR1 U2AF1-RS1 U2AF1L1 U2AF1P U2AF1RS1 U2AFBPL
Immunogen	Synthesized peptide derived from human U2AFL AA range: 131-181
Specificity	This antibody detects endogenous levels of human U2AFL
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	U2AFL;ZRSR1 U2AF1-RS1 U2AF1L1 U2AF1P U2AF1RS1 U2AFBPL;U2 small nuclear ribonucleoprotein auxiliary factor 35 kDa subunit-related protein 1 (CCCH type zinc finger;RNA-binding motif and serine/arginine rich protein 1) (U2(RNU2) small nuclear RNA auxiliary
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
