



PLB1 rabbit pAb

Catalog No	YP-Ab-09857
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PLB1 PLB
Immunogen	Synthesized peptide derived from human PLB1 AA range: 72-122
Specificity	This antibody detects endogenous levels of human PLB1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PLB;PLB/LIP;PLB1;PLB1 PLB;Phospholipase B1;membrane-associated (Phospholipase B) (hPLB) (Phospholipase B/lipase) (PLB/LIP) [Includes: Phospholipase A2 (EC 3.1.1.4);Lysophospholipase (EC 3.1.1.5)]
Observed Band	170kDa
Calculated Molecular Weight	163kDa
Cell Pathway	Apical cell membrane ; Single-pass type I membrane protein . Present in the intestinal brush border membranes. .
Tissue Specificity	Expressed in the epidermis (at protein level).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
