



SSPN rabbit pAb

Catalog No	YP-Ab-09866
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	SSPN Krag
Immunogen	Synthesized peptide derived from human SSPN AA range: 154-204
Specificity	This antibody detects endogenous levels of human SSPN
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DAGA5;Kirsten ras associated protein;Krag;nanospan;NSPN;Sarcospan;SPN1;SPN2;SSPN;SSPN Krag;Sarcospan (K-ras oncogene-associated protein) (Kirsten-ras-associated protein)
Observed Band	15-20 kDa, 27 kDa
Calculated Molecular Weight	243 aa, 27 kDa
Cell Pathway	Cell membrane ; Multi-pass membrane protein . Cell membrane, sarcolemma . Cell junction, synapse, postsynaptic cell membrane ; Multi-pass membrane protein . Also found in myotendinous junctions and in the postsynaptic membrane of neuromuscular junctions. .
Tissue Specificity	Isoform 1 is expressed exclusively in heart and skeletal muscle. Isoform 2 is expressed in heart, skeletal muscle, thymus, prostate, testis, ovary, small intestine, colon and spleen.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images