



DOCK6 rabbit pAb

Catalog No	YP-Ab-09867
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	DOCK6 KIAA1395
Immunogen	Synthesized peptide derived from human DOCK6 AA range: 1460-1510
Specificity	This antibody detects endogenous levels of human DOCK6
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	dedicator of cytokinesis 6;DOCK6;KIAA1395;ZIR1;DOCK6 KIAA1395;Dedicator of cytokinesis protein 6
Observed Band	230 kDa
Calculated Molecular Weight	2047 aa, 230 kDa
Cell Pathway	Cytoplasm . Cytoplasm, perinuclear region . Mainly located near the cell surface.
Tissue Specificity	Widely expressed. Expressed at low level in spleen, cerebellum, hippocampus and in substantia nigra.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

