



CHDH rabbit pAb

Catalog No	YP-Ab-09880
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CHDH
Immunogen	Synthesized peptide derived from human CHDH
Specificity	This antibody detects endogenous levels of human CHDH
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Molecular Weight	594 aa, 65 kDa
Observed Band	60 kDa
Cell Pathway	Mitochondrion inner membrane .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

