



SYAM rabbit pAb

Catalog No	YP-Ab-09881
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	AARS2 AARSL KIAA1270
Immunogen	Synthesized peptide derived from human SYAM AA range: 495-545
Specificity	This antibody detects endogenous levels of human SYAM
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SYAM;AARS2 AARSL KIAA1270;Alanine--tRNA ligase;mitochondrial (EC 6.1.1.7) (Alanyl-tRNA synthetase) (AlaRS)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Mitochondrion .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

