



TRNP1 rabbit pAb

Catalog No	YP-Ab-09888
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	TRNP1 C1orf225 TRNP
Immunogen	Synthesized peptide derived from human TRNP1 AA range: 38-88
Specificity	This antibody detects endogenous levels of human TRNP1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	240331C6;C1orf225;TMF-regulated nuclear protein 1;TNRP;TRNP;TRNP1;TRNP1 C1orf225 TRNP
Observed Band	23-25 kDa
Calculated Molecular Weight	227 aa, 23 kDa
Cell Pathway	Nucleus .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
