



KLK4 rabbit pAb

Catalog No	YP-Ab-09891
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	KLK4 EMSP1 PRSS17 PSTS
Immunogen	Synthesized peptide derived from human KLK4
Specificity	This antibody detects endogenous levels of human KLK4
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ARM1;EMSP;EMSP1;Kallikrein 4;Kallikrein like protein 1;kallikrein related peptidase 4;KLK L1;KLK4;PROSTASE;PRSS17;PSTS;Serine protease 17;KLK4 EMSP1 PRSS17 PSTS;Kallikrein-4 (EC 3.4.21.-) (Enamel matrix serine proteinase 1) (Kallikrein-like protein 1) (KLK-L1) (Prostase) (Serine protease 17)
Observed Band	30 kDa
Calculated Molecular Weight	254 aa, 27 kDa
Cell Pathway	Secreted.
Tissue Specificity	Expressed in prostate.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
